



SELF-EXPANDING TRACHEAL STENT SYSTEM Y-SHAPED TRACHEAL STENT

The Micro-Tech Y-shaped Tracheal Stent aids in the treatment of malignant neoplasms in the tracheobronchial carina to provide relief to patients with strictures. The self-expanding design allows for ease of deployment.



KEY BENEFITS

SILICONE COVERING

The y-stent comes with a thin silicone covering to maximize airway space while preventing tissue ingrowth.

LOW-PROFILE DELIVERY SYSTEM

The low profile of the delivery system allows for ventilation during positioning.


FLEXIBILITY

The y-stent comes in a variety of sizes to meet a range of patient anatomies and has repositioning sutures for adjustment post-deployment.

SPECIFICATIONS

Y-SHAPED TRACHEAL STENT STERILE

Order Number	Main Body		Left Branch		Right Branch		Catheter		Package Units
	OD (mm)	Length (mm)	OD (mm)	Length (mm)	OD (mm)	Length (mm)	OD (mm)	Working Length (mm)	
TT81670	16	40	12	30	12	15	8 mm	650	1/Box
TT81680	16	40	12	30	12	20	8 mm	650	1/Box
TT81690	18	40	12	30	12	15	8 mm	650	1/Box
TT81700	18	40	12	30	12	20	8 mm	650	1/Box
TT81710	18	40	12	35	12	15	8 mm	650	1/Box
TT81720	18	40	12	35	12	20	8 mm	650	1/Box
TT81730	18	60	12	30	12	15	8 mm	650	1/Box
TT81740	18	60	12	30	12	20	8 mm	650	1/Box
TT81750	18	60	12	35	12	15	8 mm	650	1/Box
TT81760	18	60	12	35	12	20	8 mm	650	1/Box
TT81770	20	50	14	30	14	15	8 mm	650	1/Box
TT81780	20	50	14	30	14	20	8 mm	650	1/Box
TT81790	20	50	14	35	14	15	8 mm	650	1/Box
TT81800	20	50	14	35	14	20	8 mm	650	1/Box
TT81810	20	60	14	30	14	15	8 mm	650	1/Box
TT81820	20	60	14	30	14	20	8 mm	650	1/Box
TT81830	20	60	14	35	14	15	8 mm	650	1/Box
TT81840	20	60	14	35	14	20	8 mm	650	1/Box
TT81850	22	60	14	35	14	20	8 mm	650	1/Box



CAN'T FIND IT?

Additional items may be available. Contact us if you can't find what you need.



2855 Boardwalk Drive
Ann Arbor, MI 48104
877.552.4027

MTEndoscopy.com
customer@mtendoscopy.com

MPW80000 Rev. 1