



SELF-EXPANDING TRACHEAL STENT SYSTEM **Y-SHAPED TRACHEAL STENT**

The Micro-tech Y-shaped Tracheal Stent aids in the treatment of malignant neoplasms in the Tracheobronchial carina to provide relief to patients with strictures. The self-expanding design allows for ease of deployment.



KEY BENEFITS

← **SILICONE COVERING**


The Y-Stent comes with a thin silicone covering to maximize airway space while preventing tissue ingrowth.

← **LOW-PROFILE DELIVERY SYSTEM**

The low profile of the delivery system allows for ventilation during positioning.

← **FLEXIBILITY**

The y-stent comes in a variety of sizes to meet a range of patient anatomies and has repositioning sutures for adjustment post- deployment.



SPECIFICATIONS

Y-SHAPED TRACHEAL STENT

ORDER #	MAIN BODY		LEFT BRANCH		RIGHT BRANCH		CATHETER		PACKAGE UNITS
	OD	LENGTH	OD	LENGTH	OD	LENGTH	OD	LENGTH	
TT81670	16 mm	40 mm	12 mm	30 mm	12 mm	15 mm	8 mm	650 mm	1/Box
TT81680	16 mm	40 mm	12 mm	30 mm	12 mm	20 mm	8 mm	650 mm	1/Box
TT81690	18 mm	40 mm	12 mm	30 mm	12 mm	15 mm	8 mm	650 mm	1/Box
TT81700	18 mm	40 mm	12 mm	30 mm	12 mm	20 mm	8 mm	650 mm	1/Box
TT81710	18 mm	40 mm	12 mm	35 mm	12 mm	15 mm	8 mm	650 mm	1/Box
TT81720	18 mm	40 mm	12 mm	35 mm	12 mm	20 mm	8 mm	650 mm	1/Box
TT81730	18 mm	60 mm	12 mm	30 mm	12 mm	15 mm	8 mm	650 mm	1/Box
TT81740	18 mm	60 mm	12 mm	30 mm	12 mm	20 mm	8 mm	650 mm	1/Box
TT81750	18 mm	60 mm	12 mm	35 mm	12 mm	15 mm	8 mm	650 mm	1/Box
TT81760	18 mm	60 mm	12 mm	35 mm	12 mm	20 mm	8 mm	650 mm	1/Box
TT81770	20 mm	50 mm	14 mm	30 mm	14 mm	15 mm	8 mm	650 mm	1/Box
TT81780	20 mm	50 mm	14 mm	30 mm	14 mm	20 mm	8 mm	650 mm	1/Box
TT81790	20 mm	50 mm	14 mm	35 mm	14 mm	15 mm	8 mm	650 mm	1/Box
TT81800	20 mm	50 mm	14 mm	35 mm	14 mm	20 mm	8 mm	650 mm	1/Box
TT81810	20 mm	60 mm	14 mm	30 mm	14 mm	15 mm	8 mm	650 mm	1/Box
TT81820	20 mm	60 mm	14 mm	30 mm	14 mm	20 mm	8 mm	650 mm	1/Box
TT81830	20 mm	60 mm	14 mm	35 mm	14 mm	15 mm	8 mm	650 mm	1/Box
TT81840	20 mm	60 mm	14 mm	35 mm	14 mm	20 mm	8 mm	650 mm	1/Box
TT81850	22 mm	60 mm	14 mm	35 mm	14 mm	20 mm	8 mm	650 mm	1/Box

MANUFACTURED BY



DISTRIBUTED BY

