



# HILZO<sup>®</sup> Through-The-Scope Fully Covered Esophageal Stent




- TTS 10.5 Fr Delivery Device
- Fully Covered with PTFE
- 20mm Diameter with 28mm Flares
- Lengths 6 to 15 cm



All sizes delivered with 10.5 Fr TTS delivery catheter and require a 3.8mm working channel for TTS delivery

# HILZO<sup>®</sup> Through-The-Scope Partially Covered Esophageal Stent



All sizes delivered with 10.5 Fr TTS delivery catheter and require a 3.8mm working channel for TTS delivery

Order #	Description	Length
EVP2006HCS2	HILZO Fully Covered TTS Esophageal Stent 28mm Flares/20mm Body	6 cm
EVP2008HCS2	HILZO Fully Covered TTS Esophageal Stent 28mm Flares/20mm Body	8 cm
EVP2010HCS2	HILZO Fully Covered TTS Esophageal Stent 28mm Flares/20mm Body	10 cm
EVP2012HCS2	HILZO Fully Covered TTS Esophageal Stent 28mm Flares/20mm Body	12 cm
EEVP2015HCS2	HILZO Fully Covered TTS Esophageal Stent 28mm Flares/20mm Body	15 cm

Order #	Description	Length
EVP2006HUS2	HILZO Partially Covered TTS Esophageal Stent 28mm Flares/20mm Body	6 cm
EVP2008HUS2	HILZO Partially Covered TTS Esophageal Stent 28mm Flares/20mm Body	8 cm
EVP2010HUS2	HILZO Partially Covered TTS Esophageal Stent 28mm Flares/20mm Body	10 cm
EVP2012HUS2	HILZO Partially Covered TTS Esophageal Stent 28mm Flares/20mm Body	12 cm
EVP2013HUS2	HILZO Partially Covered TTS Esophageal Stent 28mm Flares/20mm Body	13 cm
EVP2014HUS2	HILZO Partially Covered TTS Esophageal Stent 28mm Flares/20mm Body	14 cm
EVP2015HUS2	HILZO Partially Covered TTS Esophageal Stent 28mm Flares/20mm Body	15 cm

## FEATURES

- Silicone covered Dumbbell-shaped ends designed to reduce tissue in-growth. Body Covered with PTFE
- Dual retrieval string (proximal and distal end)
- Larger distal and proximal radiopaque markers to allow good fluoroscopic visibility.
- 10.5 Fr TTS Delivery System

## FEATURES

- Uncovered Dumbbell-shaped ends designed to reduce migration
- Body Covered with PTFE
- Dual retrieval string (proximal and distal end)
- Larger distal and proximal radiopaque markers to allow good fluoroscopic visibility.
- 10.5 Fr TTS Delivery System