

# HILZO<sup>®</sup> UES

## Fully Covered Esophageal Stent



All sizes delivered  
with 14Fr OTW delivery  
catheter under  
fluoroscopic guidance

Order #	Description	Length
<b>26MM FLARES</b>		
EVP1806FC2	HILZO Fully Covered UES Esophageal Stent 26mm Flares/18mm Body OTW	6 cm
EVP1807FC2	HILZO Fully Covered UES Esophageal Stent 26mm Flares/18mm Body OTW	7 cm
EVP1808FC2	HILZO Fully Covered UES Esophageal Stent 26mm Flares/18mm Body OTW	8 cm
EVP1809FC2	HILZO Fully Covered UES Esophageal Stent 26mm Flares/18mm Body OTW	9 cm
EVP1810FC2	HILZO Fully Covered UES Esophageal Stent 26mm Flares/18mm Body OTW	10 cm
EVP1812FC2	HILZO Fully Covered UES Esophageal Stent 26mm Flares/18mm Body OTW	12 cm
EVP1815FC2	HILZO Fully Covered UES Esophageal Stent 26mm Flares/18mm Body OTW	15 cm
<b>28MM FLARES</b>		
EVP2006FC2	HILZO Fully Covered UES Esophageal Stent 28mm Flares/20mm Body OTW	6 cm
EVP2007FC2	HILZO Fully Covered UES Esophageal Stent 28mm Flares/20mm Body OTW	7 cm
EVP2008FC2	HILZO Fully Covered UES Esophageal Stent 28mm Flares/20mm Body OTW	8 cm
EVP2009FC2	HILZO Fully Covered UES Esophageal Stent 28mm Flares/20mm Body OTW	9 cm
EVP2010FC2	HILZO Fully Covered UES Esophageal Stent 28mm Flares/20mm Body OTW	10 cm
EVP2012FC2	HILZO Fully Covered UES Esophageal Stent 28mm Flares/20mm Body OTW	12 cm
EVP2015FC2	HILZO Fully Covered UES Esophageal Stent 28mm Flares/20mm Body OTW	15 cm

## FEATURES

- Proximal Flare 0.5cm in length for added room in placements close to the UES
- Distal Flare 1.5cm to prevent migration
- Silicone covered flares designed to reduce tissue in-growth.
- Body Covered with PTFE
- Dual retrieval string (proximal and distal end)
- Larger distal and proximal radiopaque markers to allow good fluoroscopic visibility.
- 14 FR Over-the-Wire Delivery System